



What are you looking for in an ELISA-Kit?

- ✓ Superior quality
- ✓ „Made in Germany„
- ✓ CE-Mark
- ✓ Advanced usability
- ✓ Competitive pricing

Well, you've found it:

abia®

efficient • easy • better

- ✓ HIV & Hepatitis
- ✓ TORCH
- ✓ Tumor markers
- ✓ Thyroid
- ✓ Steroids
- ✓ Fertility & Prenatal
- ✓ Allergy

AB Diagnostic Systems GmbH

+49 30 208 987 10
 info@ab-ds.de
 www.ab-ds.de

LX500e Color Label Printer for patient wristbands and specimen labels

Primera Technology is a leading manufacturer of speciality printing equipment such as Primera's Signature series of slide and cassette printers for use in pathology, cytology and histology labs.

To further automate the identification process Primera offers the LX500e Color Label Printer. LX500e is a compact, full-colour, label printer designed for use in labs and other medical environments. It is ideal for printing (specimen) labels, name badges, wristbands and other types of narrow format output including many waterproof substrates.

LX500e can significantly increase the efficiency of your lab while helping to reduce the risk of misidentification of specimens. Full-colour printing adds significant impact to the printed output and helps to reduce errors by eliminating handwriting. With impressive 4800 dpi

print resolution text, graphics, colour alerts and logos along with high-resolution linear and 2D bar codes can be printed.

The new incorporated print engine doubles the maximum print speed from 25mm to 51mm per second. The printer's small size and compact footprint allows it to be installed in labs, on counters or many battery-operated medical carts that have limited space.

The LX500ec model features a built-in guillotine cutter that allows cutting of die-cut, reflective and continuous label stock.

For more information, please visit primera-healthcare.eu
 Phone: +49 611 92777-0
 Email: sales@primera.eu

