



DTM SP200 Slide Printer

Print directly onto standard slides with the new UV laser printer for precise and efficient sample processing!

Features and benefits:

- Two slide hoppers with 100 slides each
- Two output hoppers with a capacity of 20 slides, magnetically attached, can be taken along
- Print quality up to 2500 dpi
- Print speed approx. 4 s per slide
- As a UV laser printer, it does not require special laser slides, but it can print on them too
- H x W x D: 450 x 210 x 490 mm, dimensions including output hoppers
- Power consumption: 240 W
- Integrated PC with touch screen
- Built-in scanner
- Prints text, symbols, matrix codes and other characters
- Print is resistant to xylene, alcohol and decalcifying solutions
- Can be operated as stand-alone solution or with your LIS

DTM Medical GmbH
Mainzer Straße 131
65187 Wiesbaden
Deutschland
+49 611 92777-0
vertrieb@dtm-medical.eu
dtm-medical.eu

