



DTM TCP450 Tissue Cassette Printer

Print directly onto standard embedding cassettes with the new UV laser printer for precise and efficient sample processing!

Features and benefits:

- Six cassette holders with 75 cassettes each or alternatively 50 cassettes each with refillable holders
- Output holders with a capacity of 20 or 100 cassettes
- Alternatively, direct output of the printed cassettes
- Print quality up to 2500 dpi
- Print speed approx. 4 s per cassette
- As a UV laser printer, it does not require special laser cassettes
- H x W x D: 510 x 310 x 480 mm, dimensions including refillable holders for 50 cassettes
- Power consumption: 240 W
- Integrated PC with touch screen
- Built-in scanner
- Prints text, symbols, matrix codes and other characters
- Print is resistant to xylene, alcohol and decalcifying solutions
- Can be operated as stand-alone solution or with your LIS

DTM Medical GmbH
Mainzer Straße 131
65187 Wiesbaden
Deutschland
+49 611 92777-0
vertrieb@dtm-medical.eu
dtm-medical.eu

